

Histone H4 (L-Lactyl-Lys12) Rabbit mAb

Catalog No: #HL008

Package Size: #HL008-1 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Histone H4 (L-Lactyl-Lys12) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purification	>95%. Protein A/G purified from cell culture supernatant.
Applications	WB;IHC;IP
Species Reactivity	Human;Mouse;Rat
Conjugates	Unconjugated
Target Name	Histone H4
Other Names	FO108;H4;H4/n;H4F2;H4FN;HIST2H4;Histone H4;HIST1H4A;HIST2H4A
Accession No.	P62805
Calculated MW	11 kDa
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20Λ C

Application Details

WB 1:500-1:2000; IHC 1:100-1:300; IP 1:50-500

Note: This product is for in vitro research use only and is not intended for use in humans or animals.