

# Histone H4 ( L-Lactyl-Lys16) Rabbit mAb

Catalog No: #HL009

Package Size: #HL009-1 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                                         |
|--------------------|-------------------------------------------------------------------------|
| Product Name       | Histone H4 ( L-Lactyl-Lys16) Rabbit mAb                                 |
| Host Species       | Rabbit                                                                  |
| Clonality          | Monoclonal                                                              |
| Isotype            | IgG                                                                     |
| Purification       | >95%. Protein A/G purified from cell culture supernatant.               |
| Applications       | WB;IHC;IF;ChIP                                                          |
| Species Reactivity | Human;Mouse;Rat                                                         |
| Conjugates         | Unconjugated                                                            |
| Target Name        | Histone H4                                                              |
| Other Names        | FO108;H4;H4/n;H4F2;H4FN;HIST2H4;Histone H4;HIST1H4A;HIST2H4A            |
| Accession No.      | P62805                                                                  |
| Calculated MW      | 11 kDa                                                                  |
| Formulation        | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage            | Store at -20Λ C                                                         |

## Application Details

WB 1:500-1:2000; IHC 1:100-1:300; IF 1:100-1:300; IP/ChIP 1:20-200;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.