

Rabbit IgG Isotype Control

Catalog No: #IC001

Package Size: #IC001 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rabbit IgG Isotype Control
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB, ELISA, FCM, IHC-P, IHC-F, IF(IHC-P), IF(IHC-F), IF(ICC)PoO IP
Conjugates	Unconjugated
Concentration	1ug/ul
Formulation	10mM phosphate buffered saline (PBS), pH7.4
Storage	Aqueous buffered solution containing 1% BSA, 50% glycerol and 0.09% sodium azide. Store at -20°C for 12 months.

Application Details

WB Dilution:1:300-1000 ELISA Dilution:1:500-1000 FCM Dilution:1:20-100 IHC-P Dilution:1:200-400 IHC-F Dilution:1:100-500 IF(IHC-P) Dilution:1:50-200
IF(IHC-F) Dilution:1:50-200 IF(ICC) Dilution:1:50-200IPoO 1:50-200

Background

Isotype controls are negative controls that are designed to measure the amount of non-specific background signal caused by the primary antibody. Isotype Control have no specificity for the target protein, yet retain all of the non-specific characteristics of an antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.