



Rabbit IgG Isotype Control PE-Cy5.5 Conjugated

Catalog No: #IC001PE-Cy5.5

Package Size: #IC001PE-Cy5.5 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Rabbit IgG Isotype Control PE-Cy5.5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB, ELISA, FCM, IHC-P, IHC-F, IF(IHC-P), IF(IHC-F), IF(ICC)
Conjugates	conjugated
Concentration	1ug/ul
Storage	Aqueous buffered solution containing 1% BSA, 50% glycerol and 0.09% sodium azide. Store at -20°C for 12 months.

Application Details

- WB Dilution:1:300-1000
- ELISA Dilution:1:500-1000
- FCM Dilution:1:20-100
- IHC-P Dilution:1:200-400
- IHC-F Dilution:1:100-500
- IF(IHC-P) Dilution:1:50-200
- IF(IHC-F) Dilution:1:50-200
- IF(ICC) Dilution:1:50-200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.