

FITC mouse IgG2a, k isotype Control

Catalog No: #IC002aFITC

Package Size: #IC002aFITC-1 100T



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	FITC mouse IgG2a, k isotype Control
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a0O κ
Applications	FC
Conjugates	FITC
Formulation	ε• «ζ 0.2%BSA ε 0.1%Proclin 300 η£·ι 'Eη • η'Ο ε ²ζ'IAζ'A²o'Ο PH7.2-7.4
Storage	4 ε'l'o'Ο ι• Ωε δΩ• ε-

Application Details

FC 5μl/ζ~'

Note: This product is for in vitro research use only and is not intended for use in humans or animals.