

PE mouse IgG2b isotype Control

Catalog No: #IC002bPE



Package Size: #IC002bPE-1 25tests #IC002bPE-2 100tests

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	PE mouse IgG2b isotype Control
Host Species	Mouse
Clonality	Monoclonal
Clone No.	4A2B5.3
Isotype	Mouse IgG2b
Purification	Purified by Protein A.
Applications	FC IF
Concentration	10 $\times$ OI/Test
Storage	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA,Store at 4°C.DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.