

Recombinant Rabbit IgG Isotype Control

Catalog No: #IC004



Package Size: #IC004-1 50ul #IC004-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Rabbit IgG Isotype Control
Host Species	Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Applications	IHC;ICC/IF;IP;FC
Formulation	Rabbit IgG in 10mM phosphate buffered saline , pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C for short term. Store at -20°C for long term. Avoid freeze/thaw cycle.

Application Details

IHC 1:100-1:500; ICC/IF 1:100-1:500; IP 1:50-1:100; FC 1:100-1:500

Note: This product is for in vitro research use only and is not intended for use in humans or animals.