

## Chicken IgY

Catalog No: #IG601



Package Size: #IG601-1 100mg #IG601-2 200mg #IG601-3 1g

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|               |                                       |
|---------------|---------------------------------------|
| Product Name  | Chicken IgY                           |
| Purification  | Water dilution- sodium sulfate ,β 85% |
| Conjugates    | None                                  |
| Concentration | β 20 mg/ml                            |
| Storage       | -20°C                                 |

## Product Description

Source: Normal chicken egg yolk

Appearance: Liquid

Filtration: 0.2um

Buffer :10mM Phosphate buffered salinePBSpH7.4

Preservative: None

Stability: Stable for 24 months

Precautions: This product can be kept at 2-8/ C for maximum 30 days. Avoid repeated freezing and thawing. If multiple usages are needed, it is recommended to be stored separately in smaller volumes.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.