

# Human IgG

Catalog No: #IG701



Package Size: #IG701-1 100mg #IG701-2 500mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|               |                                                               |
|---------------|---------------------------------------------------------------|
| Product Name  | Human IgG                                                     |
| Purification  | Affinity chromatography or caprylic acid-ammonium sulfate,95% |
| Concentration | mg/ml                                                         |
| Storage       | -20°C                                                         |

## Product Description

Source: Normal, healthy, unimmunized rabbit serum

Appearance: Liquid

Filtration: 0.2um

Buffer :10mM Phosphate buffered salinePBSpH7.4

Preservative: None

Stability: Stable for 24 months

Precautions: This product can be kept at 2-8°C for maximum 30 days. Avoid repeated freezing and thawing. If multiple usages are needed, it is recommended to be stored separately in smaller volumes.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.