

Goat Anti-Rabbit IgG Secondary Antibody Rhodamine (TRITC) Conjugated

Catalog No: #L30114

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #L30114-1 100ul #L30114-2 1.0 ml

Description

Product Name	Goat Anti-Rabbit IgG Secondary Antibody Rhodamine (TRITC) Conjugated
Host Species	Goat
Clonality	Polyclonal
Species Reactivity	Rb
Specificity	IgG (H+L)
Conjugates	Rhodamine (TRITC)
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 10mg/ml Bovine Serum Albumin, 0.05% Sodium Azide, 50% glycerol
Storage	Store at -20 °C/1 year

Application Details

Antibody Format: Whole IgG

most applications: 1:100 - 1:200

Published Papers

el at., Astragalus polysaccharides exerts anti-infective activity by inducing human cathelicidin antimicrobial peptide LL-37 in respiratory epithelial cells. In Phytother Res. On 2018 Aug by Zhao L, Tan S et al.. PMID: 29672953, , (2018)

[PMID:29672953](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.