

Goat Anti-Rabbit IgG Secondary Antibody AF488 Conjugated

Catalog No: #L3016

Package Size: #L3016-1 100ul #L3016-2 1.0 ml

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Goat Anti-Rabbit IgG Secondary Antibody AF488 Conjugated
Host Species	Goat
Clonality	Polyclonal
Species Reactivity	Rb
Specificity	IgG (H+L)
Conjugates	AF488
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 10mg/ml Bovine Serum Albumin, 0.05% Sodium Azide, 50% glycerol
Storage	Store at -20 °C/1 year

Application Details

Antibody Format: Whole IgG

most applications: 1:100 - 1:800

Published Papers

el at., Gonadal development and molecular analysis revealed the critical window for sex differentiation, and E2 reversibility of XY-male spotted scat, *Scatophagus argus*., , (2021)

PMID:

el at., A Novel Airway-Organoid Model Based on a Nano-Self-Assembling Peptide: Construction and Application in Adenovirus Infection Studies *In Int J Nanomedicine* On 2023 Sep 15 by Yun-E Xu , Di-Shu Ao et al.. PMID:37727653, , (2023)

PMID:37727653

Note: This product is for in vitro research use only and is not intended for use in humans or animals.