

Goat Anti-Human IgG Secondary Antibody Biotin Conjugated

Catalog No: #L3024



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #L3024-1 100ul #L3024-2 1.0 ml

Description

Product Name	Goat Anti-Human IgG Secondary Antibody Biotin Conjugated
Host Species	Goat
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	IgG (H+L)
Conjugates	Biotin (long spacer)
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 10mg/ml Bovine Serum Albumin, 0.05% Sodium Azide, 50% glycerol
Storage	Store at -20 °C/1 year

Application Details

Antibody Format:Whole IgG
Working with enzyme-conjugated streptavidin
ELISA: 1:20,000 - 1:200,000
Western Blot: 1:20,000 - 1:200,000
IHC/ICC: 1:500 - 1:5,000
Flow Cytometry: 1:200 - 1:1,000
IF1:200 - 1:1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.