

# Goat Anti-Mouse IgG Secondary Antibody Cy3 Conjugated

Catalog No: #L30311

Package Size: #L30311-1 100ul #L30311-2 1.0 ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                                                                            |
|--------------------|------------------------------------------------------------------------------------------------------------|
| Product Name       | Goat Anti-Mouse IgG Secondary Antibody Cy3 Conjugated                                                      |
| Host Species       | Goat                                                                                                       |
| Clonality          | Polyclonal                                                                                                 |
| Species Reactivity | Ms                                                                                                         |
| Specificity        | IgG (H+L)                                                                                                  |
| Conjugates         | Cyanine Cy3                                                                                                |
| Formulation        | 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 10mg/ml Bovine Serum Albumin, 0.05% Sodium Azide, 50% glycerol |
| Storage            | Store at -20 °C/1 year                                                                                     |

## Application Details

Antibody Format: Whole IgG

most applications: 1:100 - 1:800

Note: This product is for in vitro research use only and is not intended for use in humans or animals.