

Goat Anti-Mouse IgG Secondary Antibody FITC Conjugated

Catalog No: #L30313

Package Size: #L30313-1 100ul #L30313-2 1.0 ml

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Goat Anti-Mouse IgG Secondary Antibody FITC Conjugated
Host Species	Goat
Clonality	Polyclonal
Species Reactivity	Ms
Specificity	IgG (H+L)
Conjugates	Fluorescein (FITC)
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 10mg/ml Bovine Serum Albumin, 0.05% Sodium Azide, 50% glycerol
Storage	Store at -20 °C/1 year

Application Details

Antibody Format: Whole IgG

most applications: 1:100 - 1:200

Published Papers

el at., A novel antibody targeting CD24 and hepatocellular carcinoma in vivo by near-infrared fluorescence imaging. In Immunobiology on 2015 Dec by Hua He, Xiaojie Tu et al.. PMID:26255089 , , (2015)

[PMID:26255089](#)

el at., Activation of the intrinsic apoptosis pathway contributes to the induction of apoptosis in hepatocellular carcinoma cells by valproic acid. In Oncol Lett on 2015 Feb by Weihua Yang , Xia Zhao et al.. PMID:25621063 , , (2015)

[PMID:25621063](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.