

Goat Anti-Mouse IgG Secondary Antibody AP Conjugated

Catalog No: #L3033



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #L3033-1 100ul #L3033-2 1.0 ml

Description

Product Name	Goat Anti-Mouse IgG Secondary Antibody AP Conjugated
Host Species	Goat
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Species Reactivity	Mouse
Specificity	IgG (H+L)
Immunogen Description	Mouse IgG
Conjugates	AP
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide
Storage	Store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

Antibody Format:Whole IgG
ELISA: 1:5,000-1:50,000Western Blot: 1:5,000-1:50,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.