

Goat Anti-Mouse IgG Secondary Antibody Biotin Conjugated

Catalog No: #L3034

Package Size: #L3034-1 100ul #L3034-2 1.0 ml

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Goat Anti-Mouse IgG Secondary Antibody Biotin Conjugated
Host Species	Goat
Clonality	Polyclonal
Species Reactivity	Ms
Specificity	IgG (H+L)
Conjugates	Biotin (long spacer)
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 10mg/ml Bovine Serum Albumin, 0.05% Sodium Azide, 50% glycerol
Storage	Store at -20 °C/1 year

Application Details

Antibody Format: Whole IgG

Working with enzyme-conjugated streptavidin

ELISA: 1:20,000 - 1:200,000

Western Blot: 1:20,000 - 1:200,000

IHC/ICC: 1:500 - 1:5,000

Flow Cytometry: 1:200 - 1:1,000

IF1:200 - 1:1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.