

Goat Anti-Mouse IgG Secondary Antibody AF647 Conjugated

Catalog No: #L3038

Package Size: #L3038 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Goat Anti-Mouse IgG Secondary Antibody AF647 Conjugated
Clonality	Polyclonal
Species Reactivity	Ms
Specificity	IgG (H+L)
Conjugates	AF647
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 10mg/ml Bovine Serum Albumin, 0.05% Sodium Azide, 50% glycerol
Storage	Store at -20 °C/1 year

Application Details

Antibody Format:Whole IgG

most applications: 1:100 - 1:800

Published Papers

el at., MicroRNA-22 Inhibits the Apoptosis of Vascular Smooth Muscle Cell by Targeting p38MAPK α in Vascular Remodeling of Aortic Dissection. In Mol Ther Nucleic Acids on 2020 Aug 25 by Yu Xiao, Yudong Sun, et al..PMID:33294292, , (2020)

[PMID:33294292](#)

el at., A Novel Airway-Organoid Model Based on a Nano-Self-Assembling Peptide: Construction and Application in Adenovirus Infection StudiesIn Int J NanomedicineOn2023 Sep 16byYun-E Xu , Di-Shu Ao et al..PMID:37727654, , (2023)

[PMID:37727654](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.