

Rabbit anti-Human IgG secondary antibody HRP-conjugated

Catalog No: #L3052

Package Size: #L3052 1ml

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rabbit anti-Human IgG secondary antibody HRP-conjugated
Clonality	Polyclonal
Applications	ELISA WB IHC
Species Reactivity	Hu
Immunogen Description	The antiserum was developed in rabbit using the human IgG as the immunogen.
Conjugates	HRP
Target Name	Human IgG
Formulation	HRP-labeled Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 50% glycerol.
Storage	Store at -20 °C/1 year

Application Details

ELISA: 1:5,000 - 1:20,000

WB: 1:5,000 - 1:20,000

IHC: 1:500 - 1:1,000

Published Papers

Loris Riccardo Lopetuso, Valentina Petito, Valerio Cufinob el at., Locally injected Infliximab ameliorates murine DSS colitis: Differences in serum and intestinal levels of drug between healthy and colitic mice, L.R. Lopetuso et al. / Digestive and Liver Disease, 45(12):1017-21(2013)

[PMID:23911613](#)

el at., Locally injected Infliximab ameliorates murine DSS colitis: differences in serum and intestinal levels of drug between healthy and colitic mice. In Dig Liver Dis on 2013 Dec by Loris Riccardo Lopetuso, Valentina Petito, et al..PMID: 23911613, , (2013)

[PMID:23911613](#)

el at., Renal Cell Carcinoma With Synchronous Ipsilateral Urothelial Carcinoma of the Renal Pelvis.In Oncol Lett on 2017 Jun by Qun Lu , Junlong Zhuang,et al..PMID: 28588718, , (2017)

[PMID:28588718](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.