

Goat Anti-Mouse IgG (H&L) HRP

Catalog No: #L35007



Package Size: #L35007-1 100 ul #L35007-2 500 ul #L35007-3 1 ml

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Goat Anti-Mouse IgG (H&L) HRP
Host Species	Goat
Clonality	Polyclonal
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, WB, IH, IC
Species Reactivity	Mouse
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Mouse IgG. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Mouse IgG
Conjugates	HRP
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)

Published Papers

el at., Abundant mRNA m1A modification in dinoflagellates: a new layer of gene regulation. In EMBO Rep on 2024 Nov
by Chongping Li, Ying Li,et al..PMID:39223385, , (2024)
[PMID:39223385](https://pubmed.ncbi.nlm.nih.gov/39223385/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.