

Rabbit Anti-Mouse IgG (H&L) HRP

Catalog No: #L35025



Package Size: #L35025-1 100 ul #L35025-2 500 ul #L35025-3 1 ml

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rabbit Anti-Mouse IgG (H&L) HRP
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, WB, IH, IC
Species Reactivity	Mouse
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Mouse IgG. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Mouse IgG
Conjugates	HRP
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)

Published Papers

Liao Zhifang, Liu Xiao, Li Linxuan, Li Sikai, Xing Xingxing, Zheng Xiwen, Song Wenyu, Gui Pin, Liu Qi, Rong Guanghong, Shao Yiming, Zou Mingzhi, Liao Hongbo, Wu Xin et al., Mechanism of the Proprietary Chinese Medicine β -JiuLiWan β • to Treat Ulcerative Colitis Revealed by Network Pharmacology, Molecular Docking, and Experimental β β , ACS omega, (2025)

[PMID:40415848](https://pubmed.ncbi.nlm.nih.gov/40415848/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.