

# Rabbit Anti-Duck IgY (H&L) FITC

Catalog No: #L35080

Package Size: #L35080-1 100 ul #L35080-2 500 ul #L35080-3 1 ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                                                                                                         |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | Rabbit Anti-Duck IgY (H&L) FITC                                                                                                                         |
| Host Species          | Rabbit                                                                                                                                                  |
| Clonality             | Polyclonal                                                                                                                                              |
| Purification          | The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.                                        |
| Applications          | IF, FC                                                                                                                                                  |
| Species Reactivity    | Dk                                                                                                                                                      |
| Specificity           | By immunoelectrophoresis and ELISA this antibody reacts specifically with Duck IgY. No antibody was detected against non immunoglobulin serum proteins. |
| Immunogen Description | Duck IgY                                                                                                                                                |
| Conjugates            | FITC                                                                                                                                                    |
| Formulation           | Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide                                                  |

## Application Details

IF (1/250 - 1/2000), FC (1/250 - 1/2000)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.