

Goat Anti-Chicken IgY (H&L) - AF488

Catalog No: #L35183



Package Size: #L35183 100 ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Goat Anti-Chicken IgY (H&L) - AF488
Host Species	Goat
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, IH, IF, FC
Species Reactivity	C
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Chicken IgY. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Chicken IgY
Conjugates	AF488
Concentration	1 mg/ml
Formulation	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide

Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Published Papers

el at., Lactate modulates microglial inflammatory responses after oxygen-glucose deprivation through HIF-1α-mediated inhibition of NF-kbInBrain Res BullOn2023 AprbyYuanyuan Zhang?1,?Pengyu Jia et al..PMID: 36746287, , (2023)
[PMID:36746287](https://pubmed.ncbi.nlm.nih.gov/36746287/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.