

Rabbit Anti-Human IgG F(ab')₂ - AF532

Catalog No: #L35268



Package Size: #L35268-1 100 ul #L35268-2 500 ul #L35268-3 1 ml

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rabbit Anti-Human IgG F(ab') ₂ - AF532
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, IH, IF, FC
Species Reactivity	Human
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgG F(ab') ₂ . No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Human IgG F(ab') ₂
Conjugates	AF60• §ε η'
Formulation	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.