

Goat Anti-Rabbit IgG (H&L) - AF568

Catalog No: #L35322



Package Size: #L35322-1 100 ul #L35322-2 500 ul #L35322-3 1 ml

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Goat Anti-Rabbit IgG (H&L) - AF568
Host Species	Goat
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, IH, IF, FC
Species Reactivity	Rabbit
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Rabbit IgG. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Rabbit IgG
Conjugates	AF680, Alexa 488, FITC
Formulation	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Published Papers

et al., MicroRNA-22 Inhibits the Apoptosis of Vascular Smooth Muscle Cell by Targeting p38MAPKα in Vascular Remodeling of Aortic Dissection. In Mol Ther Nucleic Acids on 2020 Aug 25 by Yu Xiao, Yudong Sun, et al..PMID:33294292, , (2020)

[PMID:33294292](https://pubmed.ncbi.nlm.nih.gov/33294292/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.