

Donkey Anti-Rabbit IgG (H&L) DyLight 488

Catalog No: #L35656



Package Size: #L35656 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Donkey Anti-Rabbit IgG (H&L) DyLight 488
Host Species	Donkey
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	IF, FC
Species Reactivity	Rabbit
Immunogen Description	Rabbit IgG
Conjugates	DyLight 488
Formulation	Liquid in 0.01M Sodium phosphate, 0.25M NaCl, 50% Glycerol, pH7.2, 0.05% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for six months. Avoid freeze/thaw cycles.

Application Details

IF (1/100 - 1/800), FC (1/100 - 1/800)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.