

Rabbit Anti-Human IgM

Catalog No: #L35658



Package Size: #L35658-1 100ul #L35658-2 500ul #L35658-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Rabbit Anti-Human IgM
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, WB, IH
Species Reactivity	Human
Immunogen Description	Human IgM
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH7.2, 0.05% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.