

Goat Anti-Mouse IgG2a - HRP

Catalog No: #L35660



Package Size: #L35660 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Goat Anti-Mouse IgG2a - HRP
Host Species	Goat
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, WB, IH, IC
Species Reactivity	Mouse
Immunogen Description	Mouse IgG2a
Conjugates	HRP
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Proclin
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

E (1/5000 - 1/20000),

WB (1/5000 - 1/20000),

IH (1/500 - 1/1000),

IC (1/500 - 1/1000)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.