

Goat Anti-Mouse IgG2b - HRP

Catalog No: #L35663

Package Size: #L35663 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Goat Anti-Mouse IgG2b - HRP
Host Species	Goat
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, WB, IH, IC
Species Reactivity	Mouse
Immunogen Description	Mouse IgG2b
Conjugates	HRP
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Proclin
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

E (1/5000 - 1/20000),

WB (1/5000 - 1/20000),

IH (1/500 - 1/1000),

IC (1/500 - 1/1000)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.