

Goat Anti-Mouse IgM - HRP

Catalog No: #L35664



Package Size: #L35664-1 100ul #L35664-2 500ul #L35664-3 1ml

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Goat Anti-Mouse IgM - HRP
Host Species	Goat
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, WB, IH, IC
Species Reactivity	Mouse
Immunogen Description	Mouse IgM
Conjugates	HRP
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Proclin
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

E (1/5000 - 1/20000),

WB (1/5000 - 1/20000),

IH (1/500 - 1/1000),

IC (1/500 - 1/1000)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.