

Goat Anti-Rat IgG (H&L) - AcalephFluor594

Catalog No: #L35666



Package Size: #L35666-1 100ul #L35666-2 500ul #L35666-3 1ml

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Goat Anti-Rat IgG (H&L) - AcalephFluor594
Host Species	Goat
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	IHC,IF,FC,ELISA
Species Reactivity	Rat
Immunogen Description	Rat IgG
Conjugates	AcalephFluor594
Concentration	1 mg/ml
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

IHC (1/100 - 1/1000)

IF (1/100 - 1/1000)

FC (1/1000 - 1/4000)

ELISA (Use at an assay dependent concentration)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.