

Goat Anti Mouse IgG2a Biotin Conjugated

Catalog No: #L35667



Package Size: #L35667-1 100ul #L35667-2 500ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Goat Anti Mouse IgG2a Biotin Conjugated
Host Species	Goat
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	WB;IHC;ICC
Species Reactivity	Mouse
Immunogen Description	Mouse IgG2a
Conjugates	Biotin
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:5000 - 1:20000;IHC 1:500 - 1:1000; ICC 1:500 - 1:1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.