

Donkey Anti-Rabbit IgG (H&L) - AF647

Catalog No: #L35669

Package Size: #L35669-1 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Donkey Anti-Rabbit IgG (H&L) - AF647
Host Species	Donkey
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	ELISA;IHC;IF;FC
Species Reactivity	Rabbit
Immunogen Description	Rabbit IgG
Conjugates	AF647
Concentration	1 mg/ml
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide
Storage	Store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.