

Goat Anti-Mouse IgG1

Catalog No: #L35676



Package Size: #L35676-1 100ul #L35676-2 1ml

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Goat Anti-Mouse IgG1
Host Species	Goat
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	ELISA;WB;IHC
Species Reactivity	Mouse
Immunogen Description	Mouse IgG1
Conjugates	Unconjugated
Formulation	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH7.2, 0.05% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

ELISA: 1:5000 - 1:20000; WB: 1:5000 - 1:20000; IHC: 1:500 - 1:1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.