

Goat Anti-Mouse IgM

Catalog No: #L35678

Package Size: #L35678-1 100ul #L35678-2 500ul #L35678-3 1ml

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Goat Anti-Mouse IgM
Host Species	Goat
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	ELISA;WB;IHC
Species Reactivity	Mouse
Immunogen Description	Mouse IgM
Conjugates	Unconjugated
Formulation	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH7.2, 0.05% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

ELISA: 1:5000 - 1:20000;WB: 1:5000 - 1:20000; IHC: 1:500 - 1:1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.