

Nanoantibody anti Mouse for IP (HRP)

Catalog No: #L35679

Package Size: #L35679-1 100ul #L35679-2 1ml

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Nanoantibody anti Mouse for IP (HRP)
Host Species	Camel, VHH his-tag
Applications	WB;mlHC;IP;ELISA
Species Reactivity	Mouse
Conjugates	HRP
Target Name	IgG1/2a/2b/3 Kappa
Storage	Store at 2-8°C

Application Details

WB 1:2000-10000; IP 1:5000-20000; mlHC 200-500; ELISA 5000-50000

Images



Mouse IgG for IP was diluted at 1:3000 to detect non-reduced or reduced different forms of mouse IgG, 1ug/lane

Note: This product is for in vitro research use only and is not intended for use in humans or animals.