

Goat Anti-Human IgG F(ab')₂ - HRP

Catalog No: #L35680



Package Size: #L35680-1 100ul #L35680-2 500ul #L35680-3 1ml

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Goat Anti-Human IgG F(ab') ₂ - HRP
Host Species	Goat
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	ELISA;WB;IHC;ICC
Species Reactivity	Human
Immunogen Description	Human IgG F(ab') ₂
Conjugates	HRP
Formulation	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Proclin
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

ELISA 1:5000-20000;WB 1:5000-20000;IHC 1:500-1000;ICC 1:500-1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.