

## Goat Anti-Rat IgG (H&amp;L) - AF488

Catalog No: #L35682



Package Size: #L35682-1 100ul #L35682-2 500ul #L35682-3 1ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                                                                   |
|-----------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name          | Goat Anti-Rat IgG (H&L) - AF488                                                                                   |
| Host Species          | Goat                                                                                                              |
| Clonality             | Polyclonal                                                                                                        |
| Purification          | The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads. |
| Applications          | ELISA;IHC;IF;FC                                                                                                   |
| Species Reactivity    | Rat                                                                                                               |
| Immunogen Description | Rat IgG                                                                                                           |
| Conjugates            | AF488                                                                                                             |
| Formulation           | 1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide   |
| Storage               | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.                  |

## Application Details

IHC 1:100-1:1000;IF 1:100-1:1000;FC 1:1000-1:4000;ELISA Use at an assay dependent concentration

Note: This product is for in vitro research use only and is not intended for use in humans or animals.