

Mouse Anti Human IgG2

Catalog No: #L35684



Package Size: #L35684-1 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Mouse Anti Human IgG2
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1,Kappa
Purification	The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human
Conjugates	Unconjugated
Concentration	2mg/mL
Formulation	Liquid in PBS containing, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:1000-10000; ELISA 1:10000-20000; IHC 1:50-100; IF 1:50-1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.