

HRP*Monoclonal mouse Anti Human IgG4

Catalog No: #L35685

Package Size: #L35685-1 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	HRP*Monoclonal mouse Anti Human IgG4
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.
Applications	IF, WB, IHC, ELISA
Conjugates	Unconjugated
Concentration	2mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-15°C to -25°C/1 year (Do not lower than -25°C)

Application Details

WB 1:1000-10000,

ELISA 1:10000-20000,

IHC 1:50-100

IF 1:50-1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.