

Goat Anti-Human IgG (H&L)

Catalog No: #L35686



Package Size: #L35686-1 100ul #L35686-2 1ml

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Goat Anti-Human IgG (H&L)
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	WB;IHC;ELISA
Species Reactivity	Human
Immunogen Description	Human IgG
Formulation	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH7.2, 0.05% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

WB 1:5000-20000; IHC 1:500-1000; ELISA 1:5000-20000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.