

Goat Anti-Mouse IgG1 - AF405

Catalog No: #L35687

Package Size: #L35687 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Goat Anti-Mouse IgG1 - AF405
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	ELISA;IHC;IF;FC
Species Reactivity	Mouse
Formulation	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

IHC 1/100 - 1/1000oO IF 1/100 - 1/1000oO FC 1/1000 - 1/4000oO E Use at an assay dependent concentration

Note: This product is for in vitro research use only and is not intended for use in humans or animals.