

Donkey Anti-Goat IgG (H&L) - Cyanine3

Catalog No: #L35688

Package Size: #L35688-1 100ul #L35688-2 500ul #L35688-3 1ml

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Donkey Anti-Goat IgG (H&L) - Cyanine3
Host Species	Donkey
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	ELISA;IHC;IF;FC
Species Reactivity	Goat
Immunogen Description	Goat IgG
Conjugates	Cyanine3
Formulation	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

IHC 1:100-1000; IF 1:100-1000; FC 1:1000-4000; ELISA Use at an assay dependent concentration

Note: This product is for in vitro research use only and is not intended for use in humans or animals.