

Rabbit Anti-Hamster IgG (H&L) - HRP

Catalog No: #L35689



Package Size: #L35689-1 100ul #L35689-2 500ul #L35689-3 1ml

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Rabbit Anti-Hamster IgG (H&L) - HRP
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	ELISA;WB;IHC;ICC
Species Reactivity	Hamster
Immunogen Description	Hamster IgG
Conjugates	HRP
Formulation	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Proclin
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

ELISA 1:5000-20000; WB 1:5000-20000; IHC 1:500-1000; ICC 1:500-1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.