

Goat Anti-Mouse IgM - AF594

Catalog No: #L35690



Package Size: #L35690-1 100ul #L35690-2 1ml

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Goat Anti-Mouse IgM - AF594
Host Species	Goat
Clonality	Polyclonal
Applications	IHC;IF;FC
Species Reactivity	Mouse
Immunogen Description	Mouse IgM
Conjugates	AF594
Formulation	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

IHC 1:100-1000; IF 1:100-1000; FC 1:1000-4000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.