

Mouse Anti-Human IgE - HRP

Catalog No: #L35693

Package Size: #L35693-1 100ul #L35693-2 500ul #L35693-3 1ml

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Anti-Human IgE - HRP
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	ELISA;WB;IHC;ICC
Species Reactivity	Human
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgE. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Human IgE
Conjugates	HRP
Formulation	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Proclin
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

ELISA 1:5000-20000; WB 1:5000-20000; IHC 1:500-1000; ICC 1:500-1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.