

Podoplanin Antibody

Catalog No: #P1051

Package Size: #P1051 0.1m

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Podoplanin Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	IHC650
Applications	IHC
Species Reactivity	Human
Formulation	Tris Buffer, pH 7.3 - 7.7, with 1% BSA and
Storage	Store at 2-8C. Do not freeze.

Application Details

Recommended

working dilution range:1:100 - 1:200

Images



Podoplanin [P1051] on Testicular Cancer

Product Description

Podoplanin is a transmembrane mucoprotein specifically expressed in the endothelium of lymphatic capillaries, while remaining absent from the blood vasculature. The protein is co-localized with VEGFR3/FLT4 in normal skin and kidney. Anti-Podoplanin is useful in the identification of lymphangiomas, Kaposi sarcomas, epithelioid mesotheliomas, hemangioblastomas, seminomas, and some angiosarcomas which likely have lymphatic differentiation.

References

- Ordenez NG. Adv Anat Pathol. 2006; 13:83-8.
- Ordenez NG. Hum Pathol. 2005; 36:372-80.
- Niakosari F, et al. Arch Dermatol. 2005; 141:440-4.
- Galambos C, et al. Pediatr Dev Pathol. 2005; 8:191-9.
- Fukunaga M. Histopathology. 2005; 46:396-402.

Fogt F, et al. Oncol Rep. 2004; 11:47-50.
Kahn HJ, et al. Mod Pathol. 2002; 15:434-40.
Franke FE, et al. J Cutan Pathol. 2004; 31:362-7.
Kalof AN, et al. Adv Anat Pathol. 2009; 16:62-4.
Chu AY, et al. Mod Pathol. 2005; 18:105-10.
Iczkowski KA, et al. Hum Pathol. 2008; 39:275C81.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.