

iCRT-5

Catalog No: #S08713



Package Size: #S08713-1 1 mg #S08713-2 10 mg #S08713-3 100 mg #S08713-4 25 mg #S08713-5 5 mg #S08713-6 50 mg
 Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	iCRT-5
Other Names	β catenin;Wnt;Wnt/betacatenin;Wnt/ β -catenin;Wnt/b-catenin;iCRT5;iCRT-5;iCRT5;Apoptosis;bcatenin;betacatenin;beta-catenin
Calculated MW	367.44
Storage	Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.

Application Details

SMILES: CC(=O)N(CCC(=O)SC(=O)N)C(=O)O
 MFC: MFCD00703657
 Concentration: 30 mg/mL (81.65 mM) in DMSO
 Application: Apoptosis|Wnt/beta-catenin

Note: This product is for in vitro research use only and is not intended for use in humans or animals.