

IDD388

Catalog No: #S08714



Package Size: #S08714-1 1 mg #S08714-2 10 mg #S08714-3 100 mg #S08714-4 200 mg #S08714-5 25 mg #S08714-6 5 mg #S08714-7 50 mg
Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	IDD388
Other Names	IDD 388;IDD388;IDD-388;RNA Synthesis;RNASynthesis;DNA Synthesis;DNA/RNA Synthesis;DNASynthesis
Calculated MW	416.63
Storage	keep away from moisture . Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.

Application Details

ε ε• εO• :C16H12BrCIFNO4Smile:C(NCC1=C(F)C=C(Br)C=C1)(=O)C2=C(OCC(O)=O)C=C(Cl)C=C2η1—ε!|:99.56%ζ!Aθ\$£ε!|:DMSO:40 mg/mL (96.01 mM)|• Aη 'Hζ'»ζ §:ALR2:30 nMδΩ'ε• ·ι θ—:DNA ζ• δO ε δΩ ε • |η» θ ε ``ζ ι• Aη 'H:DNA/RNA Synthesis

Note: This product is for in vitro research use only and is not intended for use in humans or animals.