

Ilmetropium iodide

Catalog No: #S08718



Package Size: #S08718-1 1 mL #S08718-2 10 mg #S08718-3 100 mg #S08718-4 2 mg #S08718-5 200 mg #S08718-6 2.5 mg #S08718-7 5 mg #S08718-8 50 mg

Others: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	Ilmetropium iodide
Other Names	Ilmetropium iodide;Ilmetropium iodide
Calculated MW	459.36
Storage	Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.

Application Details

ε ε• εO• :C20H30INO3Smile:CCC(C1=CC=CC=C1)(CO)C(O[C@@H]2C[C@@H]([N+]3(C)C)CC[C@@H]3C2)=O.[I-]DMLε• :MFC18800419ηl—ε!}:98.0
 1%ζ1Aθ§εε1};DMSO:60 mg/mL (130.62 mM)δΩ'ε• ·i θ—:ε 'Aδ» i• 'Aη 'H:Others

Note: This product is for in vitro research use only and is not intended for use in humans or animals.