

# (S)-2-Amino-3-(4-hydroxy-3,5-diiodophenyl)propanoic acid dihydrate

Catalog No: #S08722

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #S08722-1 1 g #S08722-2 1 mL #S08722-3 5 g

## Description

Product Name	(S)-2-Amino-3-(4-hydroxy-3,5-diiodophenyl)propanoic acid dihydrate
Other Names	Endogenous Metabolite;EndogenousMetabolite;3,5-Diiodotyrosine;3,5-Diiodotyrosine Dihydrate;(S) 2 Amino 3 (4 hydroxy 3,5 diiodophenyl)propanoic acid dihydrate;(S)2Amino3(4hydroxy3,5diiodophenyl)propanoic acid dihydrate;apothyrin;Inhibitor;jodgorgon;lodogorgoic acid;inhibit;ε• η' ι ζ°ι 'E
Calculated MW	469.02
Storage	store at low temperature . Powder: -20°C for 3 years . In solvent: -80°C for 1 year . Shipping with blue ice.

## Application Details

ε ε• εO• :C9H13I2NO5Smile:O.O.NC(Cc1cc(I)c(O)c(I)c1)C(O)=ODMLε• :MFCD00150275ζ'IAθ§£ε'!;DMSO:27.5 mg/mL (58.63 mM)δΩ'ε• ι θ·—:δ»£θ°'ι• Aη 'H:Endogenous Metabolite

Note: This product is for in vitro research use only and is not intended for use in humans or animals.